



TRAVEL FOR NEGATIVE LENS  
IS +/- 4mm FROM THIS POSITION

PQ Änderungsbeschreibung 20-0409

Teilekennung : 0 (ohne Auswirkung)  
Modelländerung: NEIN  
StL-Änderung : NEIN

Spezifikation Bl.1 aktualisiert

Order number: 2421-004-889-00

Specification - WDI - Adjustable and rotatable plate for ATF6 sensor  
manufacturer: WiseDeviceInc., Canada  
type: Plate-LFOVQ-NADJ-ATF  
description: Adjustable negative lens;  
360° rotation for ATF6 sensor  
weight: 0.58kg

compatible ATF6: 2421-004-274-00;2421-004-364-00

NX	EU-D	AL-T1A	US-D	US-ML	not export controlled
				PDM Status	Freigabe
PROTECTIVE NOTE "DIN ISO 16016" TO BE OBSERVED	REV	ECC	DATE	APPROVED	GENERAL TOLERANCE OF DIMENSION, FORM, POS.
	b	20-0263	10.03.20	Hornbog	ISO 2768-mH
	c	20-0409	21.04.20	Hornbog	
				BASIC TOLERANCING PRINCIPLE	TITLE
				ISO 8015	Plate-LFOVQ-NADJ-ATF
			FIRST	DATE	NAME
			ISSUE	10.03.20	Hornbogner
			CHKD	10.03.20	Neumann
DIN A 3	ALL DIMENSIONS ARE IN MM AND INCLUDE SURFACE TREATMENT				
	DRAWING NO.				SHEET
	2421-004-889-00-0002c				2
	REPLACES				OF
					2