



Order number: 2421-004-809-00

Specification - WDI - ZAxisActuator kit:

manufacturer: WiseDeviceInc., Canada
 type: ZAA-LFOVQBS
 motion type: stepper
 travel stroke: ± 4 mm
 max. speed: 10 mm/s
 max. acceleration: 100 mm/s²
 max. resolution: 0.156µm/step
 max. load: 2kg
 weight: 1kg
 note: kit includes cables (ZAA to controller)

compatible controller: 2421-004-362-00
 (as stand alone; separate part number)

| | | | | | | | | | | | |
|--|------|---------|----------|----------|--|---------------------------|--------------------------------------|------------------------------------|-----------------------------------|--|--------------|
| NX | EU-D | | AL-T1A | | US-D | | US-ML | | not export controlled | | |
| | | | | | | | PDM Status Freigabe | | - | | |
| PROTECTIVE NOTE "DIN ISO 16016" TO BE OBSERVED | REV | ECC | DATE | APPROVED | GENERAL TOLERANCE OF DIMENSION, FORM, POS. ISO 2768-mH | SURF. TREATMT | SCALE - | | MATERIAL 975980 | | |
| | b | 20-0211 | 26.02.20 | Hornbog | | | BASIC TOLERANCING PRINCIPLE ISO 8015 | | | | |
| | | | | | | FIRST DATE NAME | | TITLE Z Actuator (ZAA-LFOVQBS) Kit | | | |
| | | | | | | ISSUE 18.02.20 Hornbogner | | | | | |
| | | | | | CHKD 18.02.20 Neumann | | proposal drawing | | | | |
| | | | | | | | | | | | |
| DIN A 3 | | | | | ALL DIMENSIONS ARE IN MM AND INCLUDE SURFACE TREATMENT | | | | DRAWING NO. 2421-004-809-00-0002b | | SHEET 2 OF 2 |
| | | | | | | | REPLACES | | | | |